



PELI
BIO THERMAL™

Crēdo™
THERMAL PACKAGING SOLUTIONS



CRĒDO™ OPERATING ROOM CONTAINERS

Blood banks, hospitals, emergency medical services, and military organisations know the value of protecting blood supplies. Safely transporting materials can mean the difference between life and death and professionals need a solution that will ensure precise thermal control with the ability to be seamlessly reinstated into inventory.

At the core of every Crēdo™ Operating Room Container is our patented TIC™ system surrounding the critical payload with walls of PCM (phase change material) to maintain required temperature range and duration.

WHY CHOOSE CRĒDO™ OPERATING ROOM CONTAINER?

- Patented TIC™ system and PCM-based
- Reusable phase change technology reduces costs and eliminates the complexity and restrictions of outdated components such as wet/dry ice and gel packs
- Designed and tested to last in harsh conditions
- Quality, durable container construction



GOLDEN HOUR™ OPERATING ROOM CONTAINER

 **DURATION**
10 Hours

 **TEMPERATURE**
SERIES 4
1°C to 6°C
SERIES 22
20°C to 24°C

 **VOLUME**
2L

CRĒDO™ OPERATING ROOM CONTAINER

 **DURATION**
16 Hours

 **TEMPERATURE**
SERIES 4 - 824
1°C to 6°C

 **VOLUME**
8L




CRÉDO™ OPERATING ROOM CONTAINERS

Series	Payload Size	Litres	Exterior Dimensions WxDxH		Interior Dimensions WxDxH		Tare Weight	
			mm	inches	mm	inches	kg	lbs
Series	Series 4 ORC	2L	330 x 241 x 216	13 x 10 x 9	152 x 127 x 108	6 x 5 x 4	4.1	9
	Series 22 ORC	2L	330 x 241 x 216	13 x 10 x 9	152 x 127 x 108	6 x 5 x 4	4.1	9
Series	Series 4 824 ORC	8L	406 x 260 x 381	16 x 10 x 15	282 x 165 x 165	11 x 7 x 7	6.4	14

Whether transporting platelets, RBCs, FFP, cord blood, as well as other blood products for surgical and EMT purposes, Crêdo™ thermal protection containers offer the safest and most reliable temperature-control transport solution.

ADDITIONAL DATA

GH02 Operating Room Container Thermal Validation Test Results
(Cool +20 Deg. C Ambient; 3 Water Payload Bags; 24 hours in the Refrigerator After Freezing)



2-Litres



8-Litres



FOR MORE INFORMATION ABOUT PELI BIOTHERMAL PRODUCTS [CLICK HERE](#)

FOR MORE INFORMATION ABOUT PELI BIOTHERMAL PRODUCTS, PLEASE CONTACT:

EMEA OFFICE - INFOEMEA@PELIBIOTHERMAL.COM T: +44 (0)1525 243 770

AMERICAS OFFICE - INFOAMERICAS@PELIBIOTHERMAL.COM T: +1(763) 412-4800 TOLL FREE US: T: +(877) 537-9800

ASIA OFFICE - INFOASIA@PELIBIOTHERMAL.COM T: +65 6681 0095

PELIBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, Crêdo™, Crêdo Cube™, Crêdo ProEnvision™, Crêdo Xtreme™, Crêdo Duracube™, Crêdo ProMed™, Sherpa Systems™, CoolGuard™ and CoolPall™ are registered trademarks of Pelican Products Inc., its subsidiaries and / or affiliates

